

Ultra-Fine™ LEEP Loops



Ultra-Fine™ LEEP Loops were the world's first single use loop electrodes designed specially for LLETZ.

Ultra-Fine™ LEEP Loops have been in continuous production for more than 27 years. They continue our tradition of providing exceptional value, reliability and quality.

- **OVERCOME THE PROBLEMS** of high cost, cross infection risks and health and safety hazards commonly associated with reprocessing reusable LEEP electrodes.
- **Ultra-Fine™ LEEP LOOPS** are designed to deliver higher cutting performance with less thermal damage to the sample.
- **2.4mm (3/32") x 10cm shaft** to fit all standard UK, US and European diathermy systems.

Rocket® LEEP Range	Size		Code
Rocket® LLETZ Diathermy Loop, for single use, packed sterile in cartons 25 units	Small	14mm	R57800
	Medium	18mm	R57801
	Large	22mm	R57802
	Square	10mm	R57803
Rocket® LLETZ Ball Electrodes, for single use, packed sterile in cartons 25 units	Standard	5mm	R57804
	Small	3mm	R57807

